

Antibodies Product Datasheet

Product Name:	RGS12 Antibody
Catalog Number:	C18311
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RGS12. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	RGS12
Protein Name:	Regulator of G-protein signaling 12
Other Name:	RGS12; Regulator of G-protein signaling 12; RGS12
Human Swiss Prot No.:	O14924
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RGS12 Polyclonal Antibody detects endogenous levels of RGS12 protein.
Storage and Stability:	-20°C/1 year