

Antibodies Product Datasheet

Product Name:	ARHGAP1 Antibody
Catalog Number:	C18364
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RHG1. AA range:189-238
Concentration:	1 mg/ml
Gene Name:	ARHGAP1
Protein Name:	Rho GTPase-activating protein 1
Other Name:	ARHGAP1; CDC42GAP; RHOGAP1; Rho GTPase-activating protein 1; CDC42 GTPase-activating protein; GTPase-activating protein rhoOGAP; Rho-related small GTPase protein activator; Rho-type GTPase-activating protein 1; p50-RhoGAP
Human Swiss Prot No.:	Q07960
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARHGAP1 Polyclonal Antibody detects endogenous levels of ARHGAP1 protein.
Storage and Stability:	-20°C/1 year