

Antibodies Product Datasheet

Product Name:	ARHGAP17 Antibody
Catalog Number:	C18367
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RHG17. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	ARHGAP17
Protein Name:	Rho GTPase-activating protein 17
Other Name:	ARHGAP17; RICH1; MSTP066; MSTP110; Rho GTPase-activating protein 17; Rho-type GTPase-activating protein 17; RhoGAP interacting with CIP4 homologs protein 1; RICH-1
Human Swiss Prot No.:	Q68EM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARHGAP17 Polyclonal Antibody detects endogenous levels of ARHGAP17 protein.
Storage and Stability:	-20°C/1 year