

Antibodies Product Datasheet

Product Name:	ARHGAP23 Antibody
Catalog Number:	C18372
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ARHGAP23 . at AA range: 470-550
Concentration:	1 mg/ml
Gene Name:	ARHGAP23
Protein Name:	Rho GTPase-activating protein 23
Other Name:	ARHGAP23; KIAA1501; Rho GTPase-activating protein 23; Rho-type GTPase-activating protein 23
Human Swiss Prot No.:	Q9P227
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARHGAP23 Polyclonal Antibody detects endogenous levels of ARHGAP23 protein.
Storage and Stability:	-20°C/1 year