

Antibodies Product Datasheet

Product Name:	ARHGEF16 Antibody
Catalog Number:	C18390
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ARHGEF16 . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	ARHGEF16
Protein Name:	Rho guanine nucleotide exchange factor 16
Other Name:	ARHGEF16; EPHEXIN4; NBR; Rho guanine nucleotide exchange factor 16; Ephexin-4
Human Swiss Prot No.:	Q5VV41
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARHGEF16 Polyclonal Antibody detects endogenous levels of ARHGEF16 protein.
Storage and Stability:	-20°C/1 year