

Antibodies Product Datasheet

Product Name:	Rhophilin Antibody
Catalog Number:	C18403
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RHPN1. AA range:474-523
Concentration:	1 mg/ml
Gene Name:	RHPN1
Protein Name:	Rhophilin-1
Other Name:	RHPN1; KIAA1929; Rhophilin-1; GTP-Rho-binding protein 1
Human Swiss Prot No.:	Q8TCX5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rhophilin Polyclonal Antibody detects endogenous levels of Rhophilin protein.
Storage and Stability:	-20°C/1 year