

Antibodies Product Datasheet

Product Name:	FAM80B Antibody
Catalog Number:	C18432
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RIMKB. AA range:177-226
Concentration:	1 mg/ml
Gene Name:	RIMKLB
Protein Name:	Beta-citryl-glutamate synthase B
Other Name:	RIMKLB; FAM80B; KIAA1238; Beta-citryl-glutamate synthase B; N-acetyl-aspartyl-glutamate synthetase B; NAAG synthetase B; NAAGS; Ribosomal protein S6 modification-like protein B
Human Swiss Prot No.:	Q9ULI2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAM80B Polyclonal Antibody detects endogenous levels of FAM80B protein.
Storage and Stability:	-20°C/1 year