

Antibodies Product Datasheet

Product Name:	RCL1 Antibody
Catalog Number:	C18444
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RCL1. AA range:324-373
Concentration:	1 mg/ml
Gene Name:	RCL1
Protein Name:	RNA 3'-terminal phosphate cyclase-like protein
Other Name:	RCL1; RNAC; RPC2; RPCL1; RTC2; HSPC338; RNA 3'-terminal phosphate cyclase-like protein
Human Swiss Prot No.:	Q9Y2P8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RCL1 Polyclonal Antibody detects endogenous levels of RCL1 protein.
Storage and Stability:	-20°C/1 year