

Antibodies Product Datasheet

Product Name:	Dvl-2 Antibody
Catalog Number:	C18494
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Dvl-2 . at AA range: 670-750
Concentration:	1 mg/ml
Gene Name:	DVL2
Protein Name:	Segment polarity protein dishevelled homolog DVL-2
Other Name:	DVL2; Segment polarity protein dishevelled homolog DVL-2; Dishevelled-2; DSH homolog 2
Human Swiss Prot No.:	O14641
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dvl-2 Polyclonal Antibody detects endogenous levels of Dvl-2 protein.
Storage and Stability:	-20°C/1 year