

Antibodies Product Datasheet

Product Name:	CD108 Antibody
Catalog Number:	C18497
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from CD108 . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	SEMA7A
Protein Name:	Semaphorin-7A
Other Name:	SEMA7A; CD108; SEMAL; Semaphorin-7A; CDw108; JMH blood group antigen; John-Milton-Hargen human blood group Ag; Semaphorin-K1; Sema K1; Semaphorin-L; Sema L; CD antigen CD108
Human Swiss Prot No.:	O75326
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD108 Polyclonal Antibody detects endogenous levels of CD108 protein.
Storage and Stability:	-20°C/1 year