

Antibodies Product Datasheet

Product Name:	Septin 8 Antibody
Catalog Number:	C18514
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SEPT8. AA range:170-219
Concentration:	1 mg/ml
Gene Name:	SEPT8
Protein Name:	Septin-8
Other Name:	SEPT8; KIAA0202; Septin-8
Human Swiss Prot No.:	Q92599
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Septin 8 Polyclonal Antibody detects endogenous levels of Septin 8 protein.
Storage and Stability:	-20°C/1 year