

Antibodies Product Datasheet

Product Name:	Nek5 Antibody
Catalog Number:	C18581
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NEK5. AA range:160-209
Concentration:	1 mg/ml
Gene Name:	NEK5
Protein Name:	Serine/threonine-protein kinase Nek5
Other Name:	NEK5; Serine/threonine-protein kinase Nek5; Never in mitosis A-related kinase 5; NimA-related protein kinase 5
Human Swiss Prot No.:	Q6P3R8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Nek5 Polyclonal Antibody detects endogenous levels of Nek5 protein.
Storage and Stability:	-20°C/1 year