

Antibodies Product Datasheet

Product Name:	Snk Antibody
Catalog Number:	C18597
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLK2. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	PLK2
Protein Name:	Serine/threonine-protein kinase PLK2
Other Name:	PLK2; SNK; Serine/threonine-protein kinase PLK2; Polo-like kinase 2; PLK-2; hPlk2; Serine/threonine-protein kinase SNK; hSNK; Serum-inducible kinase
Human Swiss Prot No.:	Q9NYY3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Snk Polyclonal Antibody detects endogenous levels of Snk protein.
Storage and Stability:	-20°C/1 year