

Antibodies Product Datasheet

Product Name:	PLK-5 Antibody
Catalog Number:	C18600
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLK5. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	PLK5
Protein Name:	Inactive serine/threonine-protein kinase PLK5
Other Name:	PLK5; PLK5P; FG060302; Inactive serine/threonine-protein kinase PLK5; Polo-like kinase 5; PLK-5
Human Swiss Prot No.:	Q496M5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PLK-5 Polyclonal Antibody detects endogenous levels of PLK-5 protein.
Storage and Stability:	-20°C/1 year