

Antibodies Product Datasheet

Product Name:	TNNI3K Antibody
Catalog Number:	C18622
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TNNI3K. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	TNNI3K
Protein Name:	Serine/threonine-protein kinase TNNI3K
Other Name:	TNNI3K; CARK; Serine/threonine-protein kinase TNNI3K; Cardiac ankyrin repeat kinase; Cardiac troponin I-interacting kinase; TNNI3-interacting kinase
Human Swiss Prot No.:	Q59H18
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TNNI3K Polyclonal Antibody detects endogenous levels of TNNI3K protein.
Storage and Stability:	-20°C/1 year