

Antibodies Product Datasheet

Product Name:	PP2A-B56- α Antibody
Catalog Number:	C18643
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP2R5A. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	PPP2R5A
Protein Name:	Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit alpha isoform
Other Name:	PPP2R5A; Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit alpha isoform; PP2A B subunit isoform B'-alpha; PP2A B subunit isoform B56-alpha; PP2A B subunit isoform PR61-alpha; PR61alpha; PP2A B subunit isoform R5-alpha
Human Swiss Prot No.:	Q15172
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PP2A-B56- α Polyclonal Antibody detects endogenous levels of PP2A-B56- α protein.
Storage and Stability:	-20°C/1 year