

Antibodies Product Datasheet

Product Name:	PR48 Antibody
Catalog Number:	C18654
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP2R3B. AA range:364-413
Concentration:	1 mg/ml
Gene Name:	PPP2R3B
Protein Name:	Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta
Other Name:	PPP2R3B; PPP2R3L; Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta; PP2A subunit B isoform PR48; Protein phosphatase 2A 48 kDa regulatory subunit
Human Swiss Prot No.:	Q9Y5P8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PR48 Polyclonal Antibody detects endogenous levels of PR48 protein.
Storage and Stability:	-20°C/1 year