

Antibodies Product Datasheet

Product Name:	Megsin Antibody
Catalog Number:	C18683
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SERPINB7. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	SERPINB7
Protein Name:	Serpin B7
Other Name:	SERPINB7; Serpin B7; Megsin; TP55
Human Swiss Prot No.:	O75635
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Megsin Polyclonal Antibody detects endogenous levels of Megsin protein.
Storage and Stability:	-20°C/1 year