

Antibodies Product Datasheet

Product Name:	SH2D5 Antibody
Catalog Number:	C18703
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SH2D5. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	SH2D5
Protein Name:	SH2 domain-containing protein 5
Other Name:	SH2D5; SH2 domain-containing protein 5
Human Swiss Prot No.:	Q6ZV89
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SH2D5 Polyclonal Antibody detects endogenous levels of SH2D5 protein.
Storage and Stability:	-20°C/1 year