

Antibodies Product Datasheet

Product Name:	APS Antibody
Catalog Number:	C18704
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human APS.
Concentration:	1 mg/ml
Gene Name:	SH2B2
Protein Name:	SH2B adapter protein 2
Other Name:	SH2B2; APS; SH2B adapter protein 2; Adapter protein with pleckstrin homology and Src homology 2 domains; SH2 and PH domain-containing adapter protein APS
Human Swiss Prot No.:	O14492
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	APS Polyclonal Antibody detects endogenous levels of APS protein.
Storage and Stability:	-20°C/1 year