

Antibodies Product Datasheet

Product Name:	Fish Antibody
Catalog Number:	C18713
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human Fish.
Concentration:	1 mg/ml
Gene Name:	SH3PXD2A
Protein Name:	SH3 and PX domain-containing protein 2A
Other Name:	SH3PXD2A; FISH; KIAA0418; SH3MD1; TKS5; SH3 and PX domain-containing protein 2A; Adapter protein TKS5; Five SH3 domain-containing protein; SH3 multiple domains protein 1; Tyrosine kinase substrate with five SH3 domains
Human Swiss Prot No.:	Q5TCZ1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fish Polyclonal Antibody detects endogenous levels of Fish protein.
Storage and Stability:	-20°C/1 year