

Antibodies Product Datasheet

Product Name:	SH3BGR Antibody
Catalog Number:	C18718
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SH3BGR. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	SH3BGR
Protein Name:	SH3 domain-binding glutamic acid-rich protein
Other Name:	SH3BGR; SH3 domain-binding glutamic acid-rich protein; SH3BGR protein; 21-glutamic acid-rich protein; 21-GARP
Human Swiss Prot No.:	P55822
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SH3BGR Polyclonal Antibody detects endogenous levels of SH3BGR protein.
Storage and Stability:	-20°C/1 year