

Antibodies Product Datasheet

Product Name:	Endophilin I Antibody
Catalog Number:	C18730
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Endophilin I . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	SH3GL2
Protein Name:	Endophilin-A1
Other Name:	SH3GL2; CNSA2; SH3D2A; Endophilin-A1; EEN-B1; Endophilin-1; SH3 domain protein 2A; SH3 domain-containing GRB2-like protein 2
Human Swiss Prot No.:	Q99962
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Endophilin I Polyclonal Antibody detects endogenous levels of Endophilin I protein.
Storage and Stability:	-20°C/1 year