

Antibodies Product Datasheet

Product Name:	Sck Antibody
Catalog Number:	C18733
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SHC2. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	SHC2
Protein Name:	SHC-transforming protein 2
Other Name:	SHC2; SCK; SHCB; SHC-transforming protein 2; Protein Sck; SHC-transforming protein B; Src homology 2 domain-containing-transforming protein C2; SH2 domain protein C2
Human Swiss Prot No.:	P98077
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Sck Polyclonal Antibody detects endogenous levels of Sck protein.
Storage and Stability:	-20°C/1 year