

Antibodies Product Datasheet

Product Name:	CD172g Antibody
Catalog Number:	C18747
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SIRPG. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	SIRPG
Protein Name:	Signal-regulatory protein gamma
Other Name:	SIRPG; SIRPB2; Signal-regulatory protein gamma; SIRP-gamma; CD172 antigen-like family member B; Signal-regulatory protein beta-2; SIRP-b2; SIRP-beta-2; CD antigen CD172g
Human Swiss Prot No.:	Q9P1W8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD172g Polyclonal Antibody detects endogenous levels of CD172g protein.
Storage and Stability:	-20°C/1 year