

Antibodies Product Datasheet

Product Name:	MIP-3 β Antibody
Catalog Number:	C18769
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MIP-3 β . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	CCL19
Protein Name:	C-C motif chemokine 19
Other Name:	CCL19; ELC; MIP3B; SCYA19; C-C motif chemokine 19; Beta-chemokine exodus-3; CK beta-11; Epstein-Barr virus-induced molecule 1 ligand chemokine; EBI1 ligand chemokine; ELC; Macrophage inflammatory protein 3 beta; MIP-3-beta; Small-inducible
Human Swiss Prot No.:	Q99731
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MIP-3 β Polyclonal Antibody detects endogenous levels of MIP-3 β protein.
Storage and Stability:	-20°C/1 year