

Antibodies Product Datasheet

Product Name:	IP-10 Antibody
Catalog Number:	C18786
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from IP-10 . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	CXCL10
Protein Name:	C-X-C motif chemokine 10
Other Name:	CXCL10; INP10; SCYB10; C-X-C motif chemokine 10; 10 kDa interferon gamma-induced protein; Gamma-IP10; IP-10; Small-inducible cytokine B10
Human Swiss Prot No.:	P02778
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IP-10 Polyclonal Antibody detects endogenous levels of IP-10 protein.
Storage and Stability:	-20°C/1 year