

Antibodies Product Datasheet

Product Name:	SLC6A8 Antibody
Catalog Number:	C18796
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC6A8. AA range:581-630
Concentration:	1 mg/ml
Gene Name:	SLC6A8
Protein Name:	Sodium- and chloride-dependent creatine transporter 1
Other Name:	SLC6A8; Sodium- and chloride-dependent creatine transporter 1; CT1; Creatine transporter 1; Solute carrier family 6 member 8
Human Swiss Prot No.:	P48029
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SLC6A8 Polyclonal Antibody detects endogenous levels of SLC6A8 protein.
Storage and Stability:	-20°C/1 year