

Antibodies Product Datasheet

Product Name:	GABA T-3 Antibody
Catalog Number:	C18799
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from GABA T-3 . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	SLC6A11
Protein Name:	Sodium- and chloride-dependent GABA transporter 3
Other Name:	SLC6A11; GABT3; GAT3; Sodium- and chloride-dependent GABA transporter 3; GAT-3; Solute carrier family 6 member 11
Human Swiss Prot No.:	P48066
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABA T-3 Polyclonal Antibody detects endogenous levels of GABA T-3 protein.
Storage and Stability:	-20°C/1 year