

Antibodies Product Datasheet

Product Name:	BTR1 Antibody
Catalog Number:	C18807
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC4A11. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	SLC4A11
Protein Name:	Sodium bicarbonate transporter-like protein 11
Other Name:	SLC4A11; BTR1; Sodium bicarbonate transporter-like protein 11; Bicarbonate transporter-related protein 1; Sodium borate cotransporter 1; NaBC1; Solute carrier family 4 member 11
Human Swiss Prot No.:	Q8NBS3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BTR1 Polyclonal Antibody detects endogenous levels of BTR1 protein.
Storage and Stability:	-20°C/1 year