

Antibodies Product Datasheet

Product Name:	Na ⁺ CP type II β Antibody
Catalog Number:	C18821
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SCN2B. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	SCN2B
Protein Name:	Sodium channel subunit beta-2
Other Name:	SCN2B; Sodium channel subunit beta-2
Human Swiss Prot No.:	O60939
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Na ⁺ CP type II β Polyclonal Antibody detects endogenous levels of Na ⁺ CP type II β protein.
Storage and Stability:	-20°C/1 year