

Antibodies Product Datasheet

Product Name:	SLC10A7 Antibody
Catalog Number:	C18826
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC10A7. AA range:168-217
Concentration:	1 mg/ml
Gene Name:	SLC10A7
Protein Name:	Sodium/bile acid cotransporter 7
Other Name:	SLC10A7; C4orf13; P7; PSEC0051; Sodium/bile acid cotransporter 7; Na(+)/bile acid cotransporter 7; Solute carrier family 10 member 7
Human Swiss Prot No.:	Q0GE19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SLC10A7 Polyclonal Antibody detects endogenous levels of SLC10A7 protein.
Storage and Stability:	-20°C/1 year