

Antibodies Product Datasheet

Product Name:	CNT2 Antibody
Catalog Number:	C18845
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC28A2. AA range:371-420
Concentration:	1 mg/ml
Gene Name:	SLC28A2
Protein Name:	Sodium/nucleoside cotransporter 2
Other Name:	SLC28A2; CNT2; Sodium/nucleoside cotransporter 2; Concentrative nucleoside transporter 2; CNT 2; hCNT2; Na(+)/nucleoside cotransporter 2; Sodium-coupled nucleoside transporter 2; Sodium/purine nucleoside co-transporter; SPNT; Solute carrier
Human Swiss Prot No.:	O43868
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CNT2 Polyclonal Antibody detects endogenous levels of CNT2 protein.
Storage and Stability:	-20°C/1 year