

Antibodies Product Datasheet

Product Name:	SLC17A2 Antibody
Catalog Number:	C18869
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC17A2. AA range:213-262
Concentration:	1 mg/ml
Gene Name:	SLC17A2
Protein Name:	Sodium-dependent phosphate transport protein 3
Other Name:	SLC17A2; NPT3; Sodium-dependent phosphate transport protein 3; Na(+)/PI cotransporter 3; Sodium/phosphate cotransporter 3; Solute carrier family 17 member 2
Human Swiss Prot No.:	O00624
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SLC17A2 Polyclonal Antibody detects endogenous levels of SLC17A2 protein.
Storage and Stability:	-20°C/1 year