

Antibodies Product Datasheet

Product Name:	SF3b130 Antibody
Catalog Number:	C18921
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SF3B3. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	SF3B3
Protein Name:	Splicing factor 3B subunit 3
Other Name:	SF3B3; KIAA0017; SAP130; Splicing factor 3B subunit 3; Pre-mRNA-splicing factor SF3b 130 kDa subunit; SF3b130; STAF130; Spliceosome-associated protein 130; SAP 130
Human Swiss Prot No.:	Q15393
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SF3b130 Polyclonal Antibody detects endogenous levels of SF3b130 protein.
Storage and Stability:	-20°C/1 year