

Antibodies Product Datasheet

Product Name:	SFRS11 Antibody
Catalog Number:	C18926
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SFRS11. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	SRSF11
Protein Name:	Serine/arginine-rich splicing factor 11
Other Name:	SRSF11; SFRS11; Serine/arginine-rich splicing factor 11; Arginine-rich 54 kDa nuclear protein; p54; Splicing factor; arginine/serine-rich 11
Human Swiss Prot No.:	Q05519
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SFRS11 Polyclonal Antibody detects endogenous levels of SFRS11 protein.
Storage and Stability:	-20°C/1 year