

Antibodies Product Datasheet

Product Name:	CD6 Antibody
Catalog Number:	C19060
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD6. AA range:250-299
Concentration:	1 mg/ml
Gene Name:	CD6
Protein Name:	T-cell differentiation antigen CD6
Other Name:	CD6; T-cell differentiation antigen CD6; T12; TP120; CD antigen CD6
Human Swiss Prot No.:	P30203
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD6 Polyclonal Antibody detects endogenous levels of CD6 protein.
Storage and Stability:	-20°C/1 year