

Antibodies Product Datasheet

Product Name:	Thrombospondin 4 Antibody
Catalog Number:	C19108
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human THBS4. AA range:561-610
Concentration:	1 mg/ml
Gene Name:	THBS4
Protein Name:	Thrombospondin-4
Other Name:	THBS4; TSP4; Thrombospondin-4
Human Swiss Prot No.:	P35443
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Thrombospondin 4 Polyclonal Antibody detects endogenous levels of Thrombospondin 4 protein.
Storage and Stability:	-20°C/1 year