

Antibodies Product Datasheet

Product Name:	BET5 Antibody
Catalog Number:	C19143
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TRAPPC1. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	TRAPPC1
Protein Name:	Trafficking protein particle complex subunit 1
Other Name:	TRAPPC1; BET5; MUM2; Trafficking protein particle complex subunit 1; BET5 homolog; Multiple myeloma protein 2; MUM-2
Human Swiss Prot No.:	Q9Y5R8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BET5 Polyclonal Antibody detects endogenous levels of BET5 protein.
Storage and Stability:	-20°C/1 year