

Antibodies Product Datasheet

Product Name:	eIF2By Antibody
Catalog Number:	C19276
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from eIF2By . at AA range: 240-320
Concentration:	1 mg/ml
Gene Name:	EIF2B3
Protein Name:	Translation initiation factor eIF-2B subunit gamma
Other Name:	EIF2B3; Translation initiation factor eIF-2B subunit gamma; eIF-2B GDP-GTP exchange factor subunit gamma
Human Swiss Prot No.:	Q9NR50
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF2By Polyclonal Antibody detects endogenous levels of eIF2By protein.
Storage and Stability:	-20°C/1 year