

Antibodies Product Datasheet

Product Name:	MPP10 Antibody
Catalog Number:	C19345
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MPP10 . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	MPHOSPH10
Protein Name:	U3 small nucleolar ribonucleoprotein protein MPP10
Other Name:	MPHOSPH10; MPP10; U3 small nucleolar ribonucleoprotein protein MPP10; M phase phosphoprotein 10
Human Swiss Prot No.:	O00566
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MPP10 Polyclonal Antibody detects endogenous levels of MPP10 protein.
Storage and Stability:	-20°C/1 year