

Antibodies Product Datasheet

Product Name:	USP48 Antibody
Catalog Number:	C19392
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human USP48. AA range:631-680
Concentration:	1 mg/ml
Gene Name:	USP48
Protein Name:	Ubiquitin carboxyl-terminal hydrolase 48
Other Name:	USP48; USP31; Ubiquitin carboxyl-terminal hydrolase 48; Deubiquitinating enzyme 48; Ubiquitin thioesterase 48; Ubiquitin-specific-processing protease 48
Human Swiss Prot No.:	Q86UV5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	USP48 Polyclonal Antibody detects endogenous levels of USP48 protein.
Storage and Stability:	-20°C/1 year