

Antibodies Product Datasheet

Product Name:	FMIP Antibody
Catalog Number:	C19460
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human THOC5. AA range:561-610
Concentration:	1 mg/ml
Gene Name:	THOC5
Protein Name:	THO complex subunit 5 homolog
Other Name:	THOC5; C22orf19; KIAA0983; THO complex subunit 5 homolog; Functional spliceosome-associated protein 79; fSAP79; NF2/meningioma region protein pK1.3; Placental protein 39.2; PP39.2; hTREX90
Human Swiss Prot No.:	Q13769
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FMIP Polyclonal Antibody detects endogenous levels of FMIP protein.
Storage and Stability:	-20°C/1 year