

Antibodies Product Datasheet

Product Name:	Versican Antibody
Catalog Number:	C19494
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VCAN. AA range:532-581
Concentration:	1 mg/ml
Gene Name:	VCAN
Protein Name:	Versican core protein
Other Name:	VCAN; CSPG2; Versican core protein; Chondroitin sulfate proteoglycan core protein 2; Chondroitin sulfate proteoglycan 2; Glial hyaluronate-binding protein; GHAP; Large fibroblast proteoglycan; PG-M
Human Swiss Prot No.:	P13611
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Versican Polyclonal Antibody detects endogenous levels of Versican protein.
Storage and Stability:	-20°C/1 year