

Antibodies Product Datasheet

Product Name:	SFT2D2 Antibody
Catalog Number:	C19499
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SFT2B. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	SFT2D2
Protein Name:	Vesicle transport protein SFT2B
Other Name:	SFT2D2; Vesicle transport protein SFT2B; SFT2 domain-containing protein 2
Human Swiss Prot No.:	O95562
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SFT2D2 Polyclonal Antibody detects endogenous levels of SFT2D2 protein.
Storage and Stability:	-20°C/1 year