

## Antibodies Product Datasheet

|                               |                                                                                                                             |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | KV $\beta$ .3 Antibody                                                                                                      |
| <b>Catalog Number:</b>        | C19548                                                                                                                      |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                           |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                      |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.<br>ELISA: 1/10000. Not yet tested in other applications. |
| <b>Modification:</b>          | Unmodified/Total                                                                                                            |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.       |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human KCNAB3. AA range:293-342                          |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                     |
| <b>Gene Name:</b>             | KCNAB3                                                                                                                      |
| <b>Protein Name:</b>          | Voltage-gated potassium channel subunit beta-3                                                                              |
| <b>Other Name:</b>            | KCNAB3; KCNA3B; Voltage-gated potassium channel subunit beta-3;<br>K(+) channel subunit beta-3; Kv-beta-3                   |
| <b>Human Swiss Prot No.:</b>  | O43448                                                                                                                      |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                     |
| <b>Specificity:</b>           | KV $\beta$ .3 Polyclonal Antibody detects endogenous levels of KV $\beta$ .3 protein.                                       |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                |