

Antibodies Product Datasheet

Product Name:	C1orf181 Antibody
Catalog Number:	C19639
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CA181. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	ZNHIT6
Protein Name:	Box C/D snoRNA protein 1
Other Name:	ZNHIT6; BCD1; C1orf181; Box C/D snoRNA protein 1; Serologically defined breast cancer antigen NY-BR-75; Zinc finger HIT domain-containing protein 6
Human Swiss Prot No.:	Q9NWK9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C1orf181 Polyclonal Antibody detects endogenous levels of C1orf181 protein.
Storage and Stability:	-20°C/1 year