

Antibodies Product Datasheet

Product Name:	CENP-M Antibody
Catalog Number:	C20447
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CENP-M
Concentration:	1 mg/ml
Gene Name:	CENPM C22orf18 ICEN39 PANE1
Protein Name:	CENP-M
Other Name:	Centromere protein M (CENP-M) (Interphase centromere complex protein 39) (Proliferation-associated nuclear element protein 1)
Human Swiss Prot No.:	Q9NSP4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CENP-M
Storage and Stability:	-20°C/1 year