

Antibodies Product Datasheet

Product Name:	A Cyclase IX Antibody
Catalog Number:	C20506
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADCY9. AA range:137-186
Concentration:	1 mg/ml
Gene Name:	ADCY9
Protein Name:	Adenylate cyclase type 9
Other Name:	ADCY9; KIAA0520; Adenylate cyclase type 9; ATP pyrophosphate-lyase 9; Adenylate cyclase type IX; Adenylyl cyclase 9
Human Swiss Prot No.:	O60503
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A Cyclase IX Polyclonal Antibody detects endogenous levels of A Cyclase IX protein.
Storage and Stability:	-20°C/1 year