

Antibodies Product Datasheet

Product Name:	Helios Antibody
Catalog Number:	C20569
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IKZF2. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	IKZF2
Protein Name:	Zinc finger protein Helios
Other Name:	IKZF2; HELIOS; ZNFN1A2; Zinc finger protein Helios; Ikaros family zinc finger protein 2
Human Swiss Prot No.:	Q9UKS7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Helios Polyclonal Antibody detects endogenous levels of Helios protein.
Storage and Stability:	-20°C/1 year